

Title	Method for molecular weight estimation of non-denatured protein in polyacrylamide gradient gel electrophoresis
Sub Title	
Author	西沢, 秀幸(Nishizawa, Hideyuki) 北, 奈津美(Kita, Natsumi) 沖村, 早希子(Okimura, Sakiko) 高尾, 絵美(Takao, Emi) 阿部, 芳廣(Abe, Yoshihiro)
Publisher	共立薬科大学
Publication year	1988
Jtitle	共立薬科大学研究年報 (The annual report of the Kyoritsu College of Pharmacy). No.33 (1988.) ,p.149- 149
JaLC DOI	
Abstract	
Notes	抄録
Genre	Technical Report
URL	https://koara.lib.keio.ac.jp/xoonips/modules/xoonips/detail.php?koara_id=AN00062898-00000033-0149

慶應義塾大学学術情報リポジトリ(KOARA)に掲載されているコンテンツの著作権は、それぞれの著作者、学会または出版社/発行者に帰属し、その権利は著作権法によって保護されています。引用にあたっては、著作権法を遵守してご利用ください。

The copyrights of content available on the Keio Associated Repository of Academic resources (KOARA) belong to the respective authors, academic societies, or publishers/issuers, and these rights are protected by the Japanese Copyright Act. When quoting the content, please follow the Japanese copyright act.

Method for Molecular Weight Estimation of Non-denatured Protein in Polyacrylamide Gradient Gel Electrophoresis*

Hideyuki NISHIZAWA, Natsumi KITA, Sakiko OKIMURA,
Emi TAKAO, and Yoshihiro ABE

西沢秀幸, 北奈津美, 沖村早希子, 高尾絵美, 阿部芳廣

An improved method to estimate molecular weights of non-denatured proteins in polyacrylamide gel electrophoresis is described. Plotting the logarithm of the relative mobilities of proteins versus gel concentrations (Ferguson-plot) produces lines, whose slopes are related to molecular weights. Using a 9 cm×9 cm×0.5 mm gradient gel (T4% C5%—T20% C5%), molecular weights of non-denatured proteins (molecular weight range of 1×10^4 — 20×10^4) can be estimated in a simple manipulation.

* 本報告は *Electrophoresis*, 9, 803—806 (1988) に発表。